

Dongguan Hust Tony Instruments Co.,Ltd 东 莞 华 科 东 尼 仪 器 有 限 公 司

Addr: No. 1 Tech 9 Road, High-New Tech Development Zone, Songshan Lake, Dongguan City 523808,China

Tel: 86-769-23010657 Fax: 86-769-23010667 Email: sale@tonyhk.hk

Drying Oven DHG-9030A



Description:

It's applicable in industrial and mining enterprises, research institutes, medical and health units for drying, Provided for desiccation, torrefaction, wax-melting and sterilization in mining industry, laboratories and scientific research institutes.



Dongguan Hust Tony Instruments Co.,Ltd

东莞华科东尼仪器有限公司

Addr: No. 1 Tech 9 Road, High-New Tech Development Zone, Songshan Lake, Dongguan City 523808,China

Tel: 86-769-23010657 Fax: 86-769-23010667 Email: sale@tonyhk.hk

Technical parameters:

Desktop: Model:DHG-9003

Desktop. Model.Dira-7003							
Model	9023A 9053A 9123A 92		9203A				
Voltage	220v,50HZ						
Power	810W 1090W 1850W 2350W						
Temperature control range	room temperature:+10~200 degrees						
Temperature fluctuation	±1 degree						
Tracking alarm	+10 degrees						
Dimensions of the working room(mm)	300x280x2	275 420x3	55x350	550x355x55	0 600x550x645		

Vertical:Model:DHG-9000

V CI CICCITI I TO CICI DI I CO O O O							
Model	9030A 9070A 9140A 9240A				9240A		
Voltage	220v,50HZ						
Power	700W 1100W 1600W 2200W						0W
Temperature control range	room temperature:+10~200 degrees						
Temperature fluctuation	±1 degree						
Tracking alarm	+10 degrees						
Dimensions of the working room(mm)	300x300x350		400x40	00x450	450x550x55	50	500x600x750

Vertical:Model:DHG-9006

vertical: Model: Dilig 7000							
Model	9036A	9076.	A	9146A		9246A	
Voltage	220v,50HZ						
Power	950W 1400W 2000W 2600W)W		
Temperature control range	room temperature:+10~300 degrees						
Temperature fluctuation	±1 degree						
Tracking alarm	+10 degrees						
Dimensions of the working room(mm)	300x300x350		400x40	0x450	450x550x55	0	500x600x750

Vertical:Model:DHG-9008

Model	9038A	9078A	9108A	9148A	9248A		
Voltage	220v,50HZ						
Power	1200W	1830W	2000W	2400W	3250W		
Temperature control range	room temperature:+20~400 degrees						
Temperature fluctuation	±1 degree						
Tracking alarm	+10 degrees						



Dongguan Hust Tony Instruments Co.,Ltd

东莞华科东尼仪器有限公司

Addr: No. 1 Tech 9 Road, High-New Tech Development Zone, Songshan Lake, Dongguan City 523808,China

Tel: 86-769-23010657 Fax: 86-769-23010667 Email: sale@tonyhk.hk

Dimensions of the working room(mm)	300x300x350	400x400x450	450x450x500	450x550x550	500x600x750

Vertical:Model:DHG-9423A/DHG-9923A/DHG-9623A/DHG-9426A/DHG-9626A

Model	DHG-9423A	DHG-9923A	DHG-9623A	DHG-9426A	DHG-9626A	
Voltage	220v,50HZ	380V,50HZ				
Power	2540W	5150W	5800W			
Temperature control range	room ter	mperature:+10~200 degrees room temperature:+10~300 degrees				
Temperature fluctuation	±1 degree					
Tracking alarm	+10 degrees					
Dimensions of the working	600x550x1300	900x800x1400	800x600x1300	600x550x1300	800x600x1300	
room(mm)						

Structural description:

1.DHG type blower drying oven is a series of products, which can be divided into table type and

vertical type according to the appearance, and divided into different specifications and varieties

according to the size of the working room and the temperature control range. The outer casing of the

drying box is made of high-quality steel plate surface sprayed with plastic, the studio is made of

high-quality stainless steel plate, and the shelf made of stainless steel wire in the room can be adjusted

arbitrarily, and the middle layer is filled with ultra-fine glass wool for heat insulation. The desktop box

door adopts double-layer tempered glass door or observation window, and the vertical box door has a

double-layer tempered glass observation window on the upper middle, so that the heating items in the

box can be clearly observed. The joint between the working room and the box door is equipped with a

high temperature resistant silicon rubber sealing ring to ensure the seal between the working room and

the box door. The power switch with indicator light, temperature controller and other operating



Dongguan Hust Tony Instruments Co.,Ltd

东莞华科东尼仪器有限公司

Addr: No. 1 Tech 9 Road, High-New Tech Development Zone, Songshan Lake, Dongguan City 523808,China

Tel: 86-769-23010657 Fax: 86-769-23010667 Email: sale@tonyhk.hk

elements are all concentrated on the control panel in front of the cabinet, which is beautiful and easy to operate.

2. The heating and constant temperature system in the oven is mainly composed of a special low-noise

motor for the oven equipped with a centrifugal impeller, an electric heater, a suitable air duct structure and a temperature sensor. When the drying box is powered on, the motor runs at the same time, and

the heat generated by the electric heater placed directly on the bottom or back of the box is discharged

upward or forward through the air duct, and then sucked into the fan after passing through the dry

items in the working room, so that the cycle continues to even out the temperature.

3.The damper regulator of the desktop and vertical large-volume blast drying box can adjust the air

volume in and out of the box in a small amount by opening the damper adjustment knob.

4. Microcomputer intelligent temperature controller adopts self-tuning PID technology. Compared with

traditional PID control, it has the characteristics of rapid temperature control, fast response, small

overshoot and high precision. Both the set temperature and the temperature in the box have digital

display at the same time, with an upper limit. Tracking alarm and digital timer function, the timing can

be up to 9999 minutes. Use light touch buttons to set parameters, easy to operate.

Features:

1.304 stainless steel, mirror polishing processing, easy to clean and maintain.

2.PID controller with over temperature alarm and timing function ensures precise and reliable



Dongguan Hust Tony Instruments Co.,Ltd 东 莞 华 科 东 尼 仪 器 有 限 公 司

Addr: No. 1 Tech 9 Road, High-New Tech Development Zone, Songshan Lake, Dongguan City 523808,China

Tel: 86-769-23010657 Fax: 86-769-23010667 Email: sale@tonyhk.hk

control.

- 3.Air circulation system with specific air flow channel ensures a good temperature uniformity performance.
- 4. Silicon door gasket with long lifetime, and easy to change.
- 5.A damper adjustment in the front ensures the gas convection enough in working chamber.