

Product Information

Laboratory Fume Hoods - Full Steel Structure, Fireproofing Deflector

Test: A fume hood provides safe ventilation to protect against exposure to hazardous or toxic fumes, vapors, or airborne particulate. A chemical fume hood is primarily used in laboratory and manufacturing applications to protect the user or environment outside the hood, but can also be used to protect the materials or experiment under the hood.

Ducted fume hood designs connect to a remote exhaust system to safely ventilate outside the building. **Ductless fume hood** configurations filter and recirculate purified air via a built-in blower and filtration system.



Internal Chamber size: 100cm*70cm
Total size: Length 150cm* Width 80cm*High 235cm
Full steel structure, fireproofing deflector.

内部尺寸:100cm*70cm
箱体外尺寸:长 150cm X 宽 80cm X 高 235cm
全钢结构,加装防导流板